

AlphaLISA® SureFire® Ultra™

# Human Tau Aggregate Specific Detection Kit

**Product number:** ALSU-ATAU-A500, ALSU-ATAU-A10K,  
ALSU-ATAU-A50K, ALSU-ATAU-A-HV



## Kit specificity:


This assay kit contains antibodies which recognize aggregated Tau. The protein detected by this kit corresponds to UniProt ID P10636. Tau is also known as Microtubule-associated protein Tau. These antibodies recognize Tau of human origin. Other species should be tested on a case-by case basis.

## Control lysate information:

Positive Control Lysate: Prepared from recombinant human Tau protein aggregate (Abcam ab246003) at a concentration of 100 ng/mL.

## Representative data:

Positive Control lysate was serially diluted in Lysis Buffer and assayed for Aggregated Tau using the Tau Aggregate SureFire Ultra kit.



The information provided in this document is for reference purposes only and may not be all-inclusive. Revvity, Inc., its subsidiaries, and/or affiliates (collectively, "Revvity") do not assume liability for the accuracy or completeness of the information contained herein. Users should exercise caution when handling materials as they may present unknown hazards. Revvity shall not be liable for any damages or losses resulting from handling or contact with the product, as Revvity cannot control actual methods, volumes, or conditions of use. Users are responsible for ensuring the product's suitability for their specific application. REVVITY EXPRESSLY DISCLAIMS ALL WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, REGARDLESS OF WHETHER ORAL OR WRITTEN, EXPRESS OR IMPLIED, ALLEGEDLY ARISING FROM ANY USAGE OF ANY TRADE OR ANY COURSE OF DEALING, IN CONNECTION WITH THE USE OF INFORMATION CONTAINED HEREIN OR THE PRODUCT ITSELF"

